# To Diff or Not to Diff: That is the question.

#### Jennifer Mills

Clinical Scientist: Haematology and Transfusion

1111

### **Session Aims**

#### 1. Why do we diff?

What are the benefits of manual differentials?

2. When should we diff?

What are the rules around manual diffs?

#### 3. Case Studies

Examples for discussion.

#### 4. Tips and Tricks

Things to consider when reviewing results.





# When should we diff? Part 1

| Blast Flag 1 <sup>st</sup><br>occasion            | New blasts flags should be checked, and a manual diff completed if abnormal cells ID'd.           |
|---|---|
| Blast flag-<br>known blasts.                      | Known blasts where the counts have changed significantly or the last diff was >7 days ago.        |
| Diff/Review<br>Flags                              | Known abnormalities where the counts have changed significantly or the last diff was >7 days ago. |
| Neutropenia<br><1.0 (1 <sup>st</sup><br>occasion) | Perform manual diff if the neutropenia does not appear genuine                                    |



## When should we diff? Part 2

Neutrophils 1.0-Perform manual diff if the neutropenia does 1.5 not appear genuine. (subsequent) Any leucocytosis e.g. monocytosis >1.5, lymphocytosis >5.0 and white cells >30 may Leukocytosis need a manual diff. Significant Any patient where counts have changed may changes need a film, if there was a film <7 days ago.



# What is a "Significant Change"?



# When do we not diff?

Stable counts, <7 days since last film

CLL

Disorders such as MDS, CML or myelofibrosis may have blasts, but are often stable and low numbers, so do not need a diff.

| Leucopenia | Patient's with total white cell counts of <1.0x10 <sup>9</sup> /L do not need regular differentials unless something else changes. |
|------------|--|
|            | CLL nationt's with stable counts does not pood   |

CLL patient's with stable counts does not need a manual differential, even if their last film was >7 days ago.



### **Case Study 1: Presentation**

A 3 year old child presents to CAUQ with the following:

- Jaundice
- Haematuria

They have no historical results. An FBC is sent to the lab, and the following results were generated.



# Case Study 1: Presenting FBC

14/02/2023 14:02 Blood Request Reason : Jaundiced, haematuria. NO

| HB          | 93    | g/L              | ( | 115   | to | 135   | ) | Auth |
|-------------|-------|------------------|---|-------|----|-------|---|------|
| WBC         | 16.2  | 10 <b>*</b> 9/L  | ( | 5.0   | to | 14.5  | ) | Auth |
| PLT         | 204   | 10 <b>*</b> 9/L  | ( | 150   | to | 410   | ) | Auth |
| RBC         | 3.23  | 10 <b>*</b> 12/L | ( | 3.90  | to | 5.30  | ) | Auth |
| НСТ         | 0.269 | L/L              | ( | 0.340 | to | 0.400 | ) | Auth |
| MCV         | 83.4  | fL               | ( | 75    | to | 87    | ) | Auth |
| MCH         | 28.8  | pg               | ( | 25.0  | to | 33.0  | ) | Auth |
| MCHC        | 346   | g/L              | ( | 315   | to | 345   | ) | Auth |
| RDW         | 12.8  |                  | ( | 11.6  | to | 14.0  | ) | Auth |
| MPV         | 9.5   | fL               | ( | 7.5   | to | 11.2  | ) | Auth |
| Neutrophils | 12.1  | 10 <b>*</b> 9/L  | ( | 1.5   | to | 8.0   | ) | Auth |
| Lymphocytes | 3.2   | 10 <b>×</b> 9/L  | ( | 2.0   | to | 8.0   | ) | Auth |
| Monocytes   | 0.7   | 10 <b>*</b> 9/L  | ( | 0.2   | to | 1.0   | ) | Auth |

#### Would you do a film/differential?



# Case Study 1: FBC Day 2

15/02/2023 09:45 Blood
Request Reason : febrile illness, haemolysis process ? cause.

| HB          | 44    | g/L              | ( 11   | ō to | 135   | ) | Auth |
|-------------|-------|------------------|--------|------|-------|---|------|
| WBC         | 13.9  | 10 <b>*</b> 9/L  | ( 5.   | ) to | 14.5  | ) | Auth |
| PLT         | 183   | 10 <b>×</b> 9/L  | ( 15   | ) to | 410   | ) | Auth |
| RBC         | 1.47  | 10 <b>*</b> 12/L | ( 3.9  | ) to | 5.30  | ) | Auth |
| НСТ         | 0.119 | L/L              | ( 0.34 | ) to | 0.400 | ) | Auth |
| MCV         | 81.0  | fL               | (7     | ō to | 87    | ) | Auth |
| MCH         | 30.0  | pg               | ( 25.  | ) to | 33.0  | ) | Auth |
| MCHC        | 371   | g/L              | ( 31   | 5 to | 345   | ) | Auth |
| RDW         | 13.0  |                  | ( 11.  | 5 to | 14.0  | ) | Auth |
| MPV         | 9.2   | fL               | (7.    | ō to | 11.2  | ) | Auth |
| Neutrophils | 7.8   | 10 <b>*</b> 9/L  | ( 1.   | ō to | 8.0   | ) | Auth |
| Lymphocytes | 4.5   | 10 <b>*</b> 9/L  | ( 2.   | ) to | 8.0   | ) | Auth |
| Monocytes   | 1.3   | 10 <b>*</b> 9/L  | ( 0.   | 2 to | 1.0   | ) | Auth |

#### Would you do a film/differential?



# Case Study 1: Repeat FBC

# A. When the WBCC >30

| D<br>  T<br>S | ate<br>ime<br>pec | 19/02/2023<br>16:28<br>HQ938414H<br>H | 19/02/2023<br>06:41<br>HQ942155G<br>H | 8 u/k<br>u/k<br>HQ942125Y<br>H | 17/02/2023<br>01:08<br>HQ935809V<br>H | 16/02/2023<br>06:30<br>HQ935668R<br>H | 15/02/2023<br>09:45<br>HQ933783D<br>H | 8 14/02/2023<br>14:02<br>HQ931200A<br>H |
|---------------|-------------------|---------------------------------------|---------------------------------------|--------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---|
|               | t                 | 70                                    | 71                                    | 90                             | C0                                    | 45                                    | 11                                    | 03                                      |
| WBC           |                   | 29.4                                  | 32.1                                  | 49.2                           | 38.9                                  | 20.4                                  | 13.9                                  | 16.2                                    |
| PLT           |                   | 118                                   | 82                                    | 91                             | 94                                    | 123                                   | 183                                   | 204                                     |
| RBC           |                   | 2.46                                  | 5.53                                  | 2.55                           | 2.08                                  | 1.40                                  | 1.47                                  | 3.23                                    |
| HCT           |                   | 0.228                                 | 0.208                                 | 0.224                          | 0.179                                 | 0.121                                 | 0.119                                 | 0.269                                   |
| MCV           |                   | 92.6                                  | 90.6                                  | 88.0                           | 85.9                                  | 86.1                                  | 81.0                                  | 83.4                                    |
| MCH           |                   | 32.0                                  | 31.0                                  | 31.3                           | 32.9                                  | 31.9                                  | 30.0                                  | 28.8                                    |
| MCH           | С.                | 345                                   | 342                                   | 356                            | 383                                   | 370                                   | 371                                   | 346                                     |
| RDW           |                   | 15.3                                  | 15.6                                  | 15.7                           | 16.0                                  | 17.8                                  | 13.0                                  | 12.8                                    |
| MPV           |                   | 9.6                                   | 10.3                                  | 9.7                            | 9.1                                   | 8.6                                   | 9.2                                   | 9.5                                     |
| NEU           |                   | 17.6                                  | 21.5                                  | ^34.4                          | ^23.9                                 | 12.0                                  | 7.8                                   | 12.1                                    |
| LYM           |                   | 10.0                                  | 8.5                                   | ^11.8                          | ^12.1                                 | 6.4                                   | 4.5                                   | 3.2                                     |

#### B. When the Lymph count >5.0

C. >7 days after the 1<sup>st</sup> film

D. You wouldn't.

#### When would you do another film/differential?

# **Case Study 1: Conclusion**

The patient has increased myeloid precursors including occasional blasts, and some atypical lymphocytes, likely associated with severe infection.

The patient's haemolysis is the most significant feature, however, a high white cell count could represent the disorder driving the haemolysis.



### **Case Study 2: Presentation**

A 31 year old male, with known Acute Myeloid Leukaemia (AML), presents to the haematology outpatient department for preclinic bloods.

An FBC is sent, and the following results are generated.



# Case Study 2: Presenting FBC

| 10/08/2022 08:20 | Blood |       |                     |   |       |    |       |   |      |
|------------------|-------|-------|---------------------|---|-------|----|-------|---|------|
| HB               |       | 92    | g/L                 | ( | 130   | to | 170   | ) | Auth |
| WBC              |       | 19.7  | 10 <b>*</b> 9/L     | ( | 4.0   | to | 11.0  | ) | Auth |
| PLT              |       | 17    | 10 <del>×</del> 9/L | ( | 150   | to | 410   | ) | Auth |
| RBC              |       | 3.08  | 10 <b>*</b> 12/L    | ( | 4.50  | to | 5.50  | ) | Auth |
| HCT              |       | 0.273 | L/L                 | ( | 0.400 | to | 0.500 | ) | Auth |
| MCV              |       | 88.8  | fL                  | ( | 83    | to | 101   | ) | Auth |
| MCH              |       | 30.0  | pg                  | ( | 27.0  | to | 32.0  | ) | Auth |
| MCHC             |       | 338   | g/L                 | ( | 315   | to | 345   | ) | Auth |
| RDW              |       | 15.5  |                     | ( | 11.6  | to | 14.0  | ) | Auth |
| MPV              |       | 7.3   | fL                  | ( | 7.5   | to | 11.2  | ) | Auth |
| Neutrophils      |       | 0.8   | 10 <del>×</del> 9/L | ( | 2.0   | to | 7.0   | ) | Auth |
| Lymphocytes      |       | 4.2   | 10 <b>*</b> 9/L     | ( | 1.0   | to | 3.0   | ) | Auth |
| Monocytes        |       | 14.6  | 10 <del>×</del> 9/L | ( | 0.2   | to | 1.0   | ) | Auth |
| Eosinophils      |       | 0.1   | 10 <b>*</b> 9/L     | ( | 0.00  | to | 0.5   | ) | Auth |
| BAS              |       | 0.0   | 10 <del>×</del> 9/L | ( | 0.0   | to | 0.1   | ) | Auth |

#### Would you do a film/differential?



# **Case Study 2: FBC History**

| Date<br>Time<br>Spec | 2 10/08/2023<br>08:20<br>HQ911687M<br>H | 2 02/08/2022<br>11:17<br>HQ908127F<br>H | 27/07/2022<br>13:48<br>HQ794291B<br>H | 21/07/2022<br>08:28<br>HQ774564C<br>H | 2 18/07/2022<br>08:46<br>HQ774572W<br>H | 01/06/2022<br>11:15<br>HQ950763K<br>H |
|----------------------|---|---|---------------------------------------|---------------------------------------|---|---------------------------------------|
| HGB                  | 92                                      | 70                                      | 74                                    | 88                                    | 93                                      | 89                                    |
| WRC                  | 19.7                                    | 4.6                                     | 1.8                                   | 2.6                                   | 1.6                                     | 1.7                                   |
| PLT                  | 17                                      | 16                                      | 39                                    | 10                                    | 14                                      | 32                                    |
| RBC                  | 3.08                                    | 2.63                                    | 2.49                                  | 2.94                                  | 3.13                                    | 2.49                                  |
| HCT                  | 0.273                                   | 0.232                                   | 0.219                                 | 0.258                                 | 0.273                                   | 0.258                                 |
| MCV                  | 88.8                                    | 88.4                                    | 87.9                                  | 87.7                                  | 87.4                                    | 103.5                                 |
| MCH                  | 30.0                                    | 30.2                                    | 29.7                                  | 29.8                                  | 29.7                                    | 36.0                                  |
| MCHC.                | 338                                     | 342                                     | 338                                   | 339                                   | 339                                     | 347                                   |
| RDW                  | 15.5                                    | 15.6                                    | 14.9                                  | 14.7                                  | 14.5                                    | 18.5                                  |
| MPV                  | 7.3                                     | 8.3                                     | 8.0                                   | 8.1                                   | 7.7                                     | 9.4                                   |
| NEU                  | 0.8                                     | 0.7                                     | 0.8                                   | 0.9                                   | 0.6                                     | 0.1                                   |
| LYM                  | 4.2                                     | 1.5                                     | 0.6                                   | 1.5                                   | 0.9                                     | 1.5                                   |
| MON                  | 14.6                                    | 2.3                                     | 0.4                                   | 0.2                                   | 0.1                                     | 0.0                                   |

#### Does this change your opinion?



# **Case Study 2: Repeat Bloods**

A film was reviewed, and a manual differential was not done. The patient was seen again 5 days later with the following results:

| 15/08/2022 11:11 Blood<br>Request Reason : AML. |                  |                               |  |
|---|------------------|-------------------------------|--|
| HB<br>WBC<br>PLT                                | 79<br>27.7<br>15 | g/L (<br>10*9/L (<br>10*9/L ( | 130 to 170 ) Auth<br>4.0 to 11.0 ) Auth<br>150 to 410 ) Auth |
| Wha   | at wo            | uld you do?                   |  |
| MCHC  | 333              | g/L (                         | 315 to 345 ) Auth  |
| RDW   | 15.7             |                               | 11.6 to 14.0 ) Auth  |
| MPV   | 9.9              | fL (                          | 7.5 to 11.2 )Auth  |
| Neutrophils                                     | ^0.4             | 10 <del>×</del> 9/L (         | 2.0 to 7.0 ) Auth  |
| Lymphocytes                                     | 2.2              | 10 <del>×</del> 9/L           | 1.0 to 3.0 ) Auth  |
| Monocytes                                       | ^24.1            | 10 <del>×</del> 9/L (         | 0.2 to 1.0 ) Auth  |

A. No film as previous film was <7 days ago

B. No film, phone and HQ.

#### C. Repeat film/diff



# **Case Study 2: Conclusion**

The patient has a blast count of 25.2x10<sup>9</sup>/L and a monocyte count of 0.0x10<sup>9</sup>/L. This represents a relapse of his known disease.

It is important to do a manual differential on known leukaemia patients with sudden increase in White cells, as this may represent a relapse.



# Case Study 3

A 61 year old female presents to her GP with the following clinical symptoms:

- Fatigue
- Easy Bruising

She has no known comorbidities, and does not drink or smoke. A FBC sample is sent to the lab.

# Case Study 3: Presenting FBC

28/01/2021 08:45 Blood Request Reason : bruising ++. NONE

| HB          | 115   | g/L              | ( 120   | to | 150   | ) | Auth |
|-------------|-------|------------------|---------|----|-------|---|------|
| WBC         | 7.4   | 10 <b>*</b> 9/L  | ( 4.0   | to | 11.0  | ) | Auth |
| PLT         | 4     | 10 <b>*</b> 9/L  | ( 150   | to | 410   | ) | Auth |
| RBC         | 3.67  | 10 <b>*</b> 12/L | ( 3.80  | to | 4.80  | ) | Auth |
| НСТ         | 0.333 | L/L              | ( 0.360 | to | 0.460 | ) | Auth |
| MCV         | 90.6  | fL               | ( 83    | to | 101   | ) | Auth |
| MCH         | 31.3  | pg               | ( 27.0  | to | 32.0  | ) | Auth |
| MCHC        | 346   | g/L              | ( 315   | to | 345   | ) | Auth |
| RDW         | 15.3  |                  | ( 11.6  | to | 14.0  | ) | Auth |
| MPV         | 10.0  | fL               | ( 7.5   | to | 11.2  | ) | Auth |
| Neutrophils | ^0.6  | 10 <b>*</b> 9/L  | ( 2.0   | to | 7.0   | ) | Auth |
| Lymphocytes | 2.8   | 10 <b>×</b> 9/L  | ( 1.0   | to | 3.0   | ) | Auth |
| Monocytes   | ^3.8  | 10 <b>*</b> 9/L  | ( 0.2   | to | 1.0   | ) | Auth |

#### Would you do a film/differential?

Flag S Comments : Variant LY MO Blast LY Blast

# Case Study 3: Presenting FBC

The patient had a manual differential, resulting in a blast count of 3.6x10<sup>9</sup>/L.

The patient had a bone marrow, genetics and flow cytometry testing, and was diagnosed with Acute Promyelocytic Leukaemia (APML)

The patient was admitted, and had daily bloods sent to the laboratory.



# Case Study 3: Repeat FBCs

| Date<br>Time<br>Spec | 31/01/2021<br>06:35<br>HQ845830P<br>H | 30/01/2021<br>18:10<br>HQ845588S<br>H | 30/01/2021<br>06:22<br>HQ845003L<br>H | 29/01/2021<br>18:23<br>HQ844630Q<br>H | 29/01/2021<br>08:30<br>HQ841967Z<br>H | 28/01/2021<br>16:44<br>HQ841462L<br>H | 28/01/202<br>08:45<br>HQ836651R<br>H |
|----------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|--------------------------------------|
|                      | on                                    | 00                                    | no                                    | 100                                   | 107                                   | 4 4 4                                 | 445                                  |
|                      | 09                                    | 90<br>11 7                            | 90<br>16 7                            | 10.0                                  |                                       |                                       |                                      |
|                      | 9.0                                   | 11.1                                  | 10.1                                  | 19.0                                  | 19.6                                  | 9.0                                   | (.4                                  |
| PLI                  | 30                                    | 22                                    | 34                                    | у                                     | 16                                    | 4                                     | 4                                    |
| RBC                  | 2.80                                  | 2.91                                  | 3.18                                  | 3.26                                  | 3.40                                  | 3.53                                  | 3.67                                 |
| НСТ                  | 0.255                                 | 0.265                                 | 0.290                                 | 0.298                                 | 0.309                                 | 0.322                                 | 0.333                                |
| MCV                  | 91.0                                  | 91.2                                  | 91.1                                  | 91.4                                  | 90.9                                  | 91.3                                  | 90.6                                 |
| MCH                  | 31.9                                  | 30.9                                  | 30.9                                  | 30.6                                  | 31.5                                  | 31.5                                  | 31.3                                 |
| MCHC.                | 350                                   | 338                                   | 339                                   | 334                                   | 347                                   | 345                                   | 346                                  |
| RDW                  | 15.3                                  | 15.4                                  | 15.4                                  | 15.0                                  | 15.0                                  | 14.9                                  | 15.3                                 |
| MPV                  | 8.5                                   | 7.9                                   | 7.9                                   | 7.8                                   | 7.5                                   | 9.4                                   | 10.0                                 |
| NEU                  |                                       |                                       |                                       |                                       | ^0.9                                  |                                       | ^0.6                                 |
| NEUT                 |                                       | 1.5                                   |                                       |                                       | 0.6                                   |                                       | 0.5                                  |
|                      |                                       |                                       |                                       |                                       |                                       |                                       |                                      |

A. No film as previous film was <7 days ago

B. Everyday

C. When the white cell count became high

#### When would you repeat the film?

### Case Study 3: More FBCs!

| Date<br>Time<br>Spec | 06/02/2021<br>06:28<br>HQ862138S<br>H | 04/02/2021<br>04:30<br>HQ856271T<br>H | 03/02/2021<br>16 <b>:</b> 10<br>HQ911264T<br>H | 03/02/2021<br>06:30<br>HQ853203N<br>H | 02/02/2021<br>06:30<br>HQ849211Y<br>H | 01/02/2021<br>20:30<br>HQ849198S<br>H | 01/02/2021<br>08:05<br>HQ846445Z<br>H |
|----------------------|---------------------------------------|---------------------------------------|--|---------------------------------------|---------------------------------------|---------------------------------------|---------------------------------------|
| HGR                  | 80                                    | 77                                    | 81   | 80                                    | 87                                    | 83                                    | 87                                    |
| WBC                  | 42.2                                  | 60.4                                  | 71.0   | 61.3                                  | 41.1                                  | 24.6                                  | 16.1                                  |
| PLT                  | 30                                    | 27                                    | 37   | 36                                    | 46                                    | 50                                    | 60                                    |
| RBC                  | 2.65                                  | 2.49                                  | 2.60   | 2.55                                  | 2.79                                  | 2.69                                  | 2.77                                  |
| НСТ                  | 0.244                                 | 0.233                                 | 0.240  | 0.235                                 | 0.258                                 | 0.247                                 | 0.251                                 |
| MCV                  | 92.0                                  | 93.5                                  | 92.4   | 92.2                                  | 92.5                                  | 91.8                                  | 90.8                                  |
| MCH                  | 30.3                                  | 31.0                                  | 31.1   | 31.2                                  | 31.2                                  | 30.9                                  | 31.3                                  |
| MCHC.                | 329                                   | 331                                   | 336  | 338                                   | 337                                   | 337                                   | 345                                   |
| RDW                  | 15.7                                  | 16.2                                  | 15.7   | 15.9                                  | 15.4                                  | 15.5                                  | 15.3                                  |
| MPV                  | 8.1                                   | 8.0                                   | 8.5  | 8.8                                   | 7.8                                   | 8.0                                   | 8.3                                   |
| NLR                  |                                       |                                       | NLRE   |                                       |                                       |                                       |                                       |
| anrb                 | ^0.0                                  | ^0.0                                  | ^0.0   | ^0.1                                  | ^0.1                                  |                                       | ^0.0                                  |
|                      |                                       |                                       |  |                                       |                                       |                                       |                                       |

#### When would you repeat the film?

A. 04/02/2021, as this is >7 days since previous film

B. Everyday

C. When the white cell count >30x10<sup>9</sup>/L

# **Case Study 3: Conclusion**

Despite treatment, this patient's WBCC varies widely. In APML this is associated with a condition called "differentiation syndrome."

In the case of APML, it is useful to do a differential every day, as this helps determine if the patient is responding to management strategies, and to ensure they are not neutropenic.

### Tips and Tricks: Deciding to Diff

| Think about the FBC.                | Think about all of the values when deciding if the FBC has changed.  |
|-------------------------------------|--|
| Know the<br>last diff.              | Has is been >7 days? What did they see on that film diff?  |
| Don't leave<br>the analyzer<br>diff | If you choose not to make a film, make sure you delete the analyser differential.                              |
| If in doubt-<br>diff!               | There's no harm in doing a diff! If it's the same as the analyser, it just means we don't have to do it again! |



# Tips and Tricks: Doing the Diff

1 blast =/= If you see one blast, one promyelocyte etc, you don't have to a diff, just remember to mention it!

Is it a blast? Don't call it atypical! If you think you see a blast, call it that! Try to avoid "atypical cells" if you can as doctors don't know what this means.

Everything is on a spectrum! Cells are in the process of development, so may not be textbook examples. Consider the features they have e.g. granules or size, and make your decision using that rationale



# Thanks for Listening

#### Any Questions?

Jennifer Mills

<u>[ennifer.mills@porthosp.nhs.uk</u>

023 92 28 5774